



Assemble Base Cell

1. PREPARE STOCK SOLUTIONS

- Add 1 mL of mineral oil to the glass vial using a 1 mL pipette.
- Add the lipids in the above table to the glass vial on top of the mineral oil using the appropriate glass syringe.
- Briefly vortex the lipid-oil mixture for 5 seconds to mix.
- Evaporate the chloroform from the lipid-oil mixture:
 - Place glass vial in a 55°C dry bath in a fume hood.
 - Shield with aluminum foil to protect from light.
 - Evaporate uncovered for 4 hours.
- In a glass bottle, add 4 mL of chloroform using a glass 10 mL serological pipette.
- Clean syringes by rinsing with chloroform 5 times into an empty glass bottle. Store the syringes with the plunger removed inside the fume hood for (3-4) hr to allow remaining chloroform to evaporate.
- After 4 hr of incubation at 55°C in a dry bath in a fume hood, allow the lipid-oil mixture to cool to room temperature for 10-15 min.
- Return the plungers to the syringes and store them in their designated location.
- Dispose of chloroform waste following applicable chemical safety guidelines.

2. ASSEMBLE OUTER SOLUTIONS

- Prepare 1.5 mL microcentrifuge tubes labeled with the appropriate reaction.
- Mix glucose stock solution and water according to the following table:

Component	Cells + deGFP DNA (μL)	Cells - deGFP DNA (μL)
Glucose (2 M)	570	570
Water	430	430
Total	1000	1000

3. ASSEMBLE CYTOSOL REACTIONS

- Remove all components listed in the Composition table above from appropriate cold storage.
- Thaw reagents on ice.
- Prepare 1.5 mL microcentrifuge tubes, on ice, to assemble reactions into.
- For a given reaction, assemble by adding the volume of reagents from the Composition table in the order listed. Pay special attention to the handling of the Cytosol components:
 - Vortex SMix: Ensure thorough mixing; 10s vortex / 10s rest on ice; should be transparent with no visible precipitate; and add to the reaction tubes.
 - Vortex or pipette mix tRNA, and add to the reaction tubes.
 - Vortex or pipette mix PMix, and add to the reaction tubes.
 - Do NOT vortex** ribosomes: *gently* pipette mix or flick the tube, and add to the reaction tubes.
 - Add remaining reactions in the order they appear
- Mix the master mix thoroughly by pipetting up and down (10-15) times until it appears homogeneous and clear.
- Close lids on the microcentrifuge tubes and briefly spin down to eliminate bubbles.

- Pipette out 10 μL of the reaction for osmolarity check using a Vapor Pressure Osmometer before starting encapsulation.
- Hold assembled reactions on ice until ready for encapsulation.

4. ENCAPSULATE CYTOSOLS INTO LIPOSOMES

- Set up a 1.5 mL tube rack with two 1.5 mL microcentrifuge tubes for each liposome encapsulation. Number the tubes according to the number of reactions assembled. Label the two tubes for each reaction:
 - T**—transfer
 - L**—liposomes
- Add 300 μL of the appropriate glucose outer solution to each of the tubes labeled **T**.
- Add 150 μL of the lipid-oil mixture (at room temperature) on top of each assembled Cytosol reaction.
- Emulsify the lipid-oil and Cytosol reaction by running the tube along a row of empty slots on the 1.5 mL tube rack. Run it down 20–30 times until the solution forms a stable emulsion with an even milky color.
- Immediately layer each emulsion over the outer solution. Slowly pipette the entire emulsion down the side of the corresponding **T** tube.
- Centrifuge **T** tubes at 9000 g for 10 min at room temperature to pellet the liposomes.
- Extract the liposomes from each **T** tube:
 - Remove the oil layer and lipid debris from the top of each **T** tube by gently pipetting with a 1000 μL pipette set to 800 μL .
 - Gently pipette mix the pellet 10-15 times with the outer solution.
 - Extract liposomes by pipetting 50-100 μL of pellet and outer solution from **T** and transfer liposome sample to the respective liposome tube **L**.
- Hold liposomes on ice until you are prepared to begin measurement.
- Pipette the liposomes into a well on a 384-well glass bottom plate. If the density appears too high under the microscope, dilute the liposomes with the outer solution for better data analysis.
- Begin measurement.
- Add a black dot to the lid of each of Cytosol component. The number of dots indicates freeze-thaw cycles.