



## Nucleus Protocols

# Make Ribosomes

## 1. PREPARE STOCK BUFFERS

- Make the following stock solutions. Use ultrapure water (18.2 M $\Omega$ , e.g., Milli-Q) and keep everything RNase-free.

Stocks	Final Concentration (mM)	MW (g/mol)	Mass to add (g)	Final Vol (mL)	Storage (°C)
HEPES-KOH (pH 7.6)	1000	238.3	238.3	1000	room temp
Ammonium Chloride	1000	53.49	53.49	1000	room temp
Potassium Chloride	1000	74.55	74.55	1000	room temp
Magnesium Acetate	1000	214.45	214.45	1000	room temp
Sodium Hydroxide	500	40.00	20	1000	room temp

- Adjust the pH of HEPES-KOH to 7.6 with Potassium Hydroxide.
- Prepare Sodium Hydroxide (0.5 M) from pellets (used to clean the HIC column).

## 2. PREPARE STABLE BUFFERS

- Ribosome Lysis Buffer** — used to resuspend the cell pellet for sonication.

Reagent	Final Concentration (mM)	Stock Concentration (mM)	Volume to Add (mL)
HEPES-KOH (pH 7.6)	10	1000	15
Magnesium Acetate	10	1000	15
Potassium Chloride	50	1000	75
TCEP	1	500	3
Ultrapure water	—	—	1392
<b>Total</b>			<b>1500</b>

- Ribosome Salting Out Buffer** — added 1:1 to clarified lysate to bring ammonium sulfate to 1.5 M and precipitate contaminants.

Reagent	Final Concentration (mM)	Stock Concentration (mM)	Volume to Add (mL)
Ammonium Sulfate	3000	n/a (dry)	198.21 g
HEPES-KOH (pH 7.6)	10	1000	5
Magnesium Acetate	10	1000	5
Potassium Chloride	50	1000	25
TCEP	1	500	1
Ultrapure water	—	—	~265.8
<b>Total</b>			<b>500</b>

- Ribosome Wash Buffer** — HIC running/wash buffer; high ammonium sulfate drives hydrophobic binding.

Reagent	Final Concentration (mM)	Stock Concentration (mM)	Volume to Add (mL)
Ammonium Sulfate	1500	n/a (dry)	396.42 g
HEPES-KOH (pH 7.6)	20	1000	40
Magnesium Acetate	10	1000	20
TCEP	1	500	4
Ultrapure water	—	—	~1539.6
<b>Total</b>			<b>2000</b>

- Ribosome Elution Buffer** — HIC elution buffer; low salt releases bound ribosomes.

Reagent	Final Concentration (mM)	Stock Concentration (mM)	Volume to Add (mL)
HEPES-KOH (pH 7.6)	20	1000	20
Magnesium Acetate	10	1000	10
TCEP	1	500	2
Ultrapure water	—	—	968
<b>Total</b>			<b>1000</b>

- Cushion Buffer** — dense sucrose cushion through which ribosomes are pelleted during ultracentrifugation.

Reagent	Final Concentration (mM)	Stock Concentration (mM)	Volume to Add (mL)
Sucrose	30% (w/v)	n/a (dry)	300 g
HEPES-KOH (pH 7.6)	20	1000	20
Ammonium Chloride	30	1000	30
Magnesium Acetate	10	1000	10
TCEP	1	500	2
Ultrapure water	—	—	~638
<b>Total</b>			<b>1000</b>

- Ribosome Buffer** — final resuspension and storage buffer for purified ribosomes.

Reagent	Final Concentration (mM)	Stock Concentration (mM)	Volume to Add (mL)
HEPES-KOH (pH 7.6)	20	1000	10
Magnesium Acetate	6	1000	3
Potassium Chloride	30	1000	15
TCEP	1	500	1
Ultrapure water	—	—	471
<b>Total</b>			<b>500</b>

### 3. PREPARE WORKING BUFFERS

- Acetic Acid (0.1 M)** — prepared by diluting glacial acetic acid; used to clean the HIC column.

Reagent	Final Concentration (mM)	Stock Concentration (mM)	Volume to Add (mL)
Acetic Acid (glacial)	100	17 400	5.75
Ultrapure water	—	—	994.25
<b>Total</b>			1000

- Ethanol (20% v/v)** — used to wash and store the HIC column.

Reagent	Final Concentration	Stock Concentration	Volume to Add (mL)
Ethanol	20% (v/v)	100% (v/v)	200
Ultrapure water	—	—	800
<b>Total</b>			<b>1000</b>

### 4. CELL CULTURE

- Add 5 mL Luria Broth (LB) under sterile conditions to two (2) 14 mL culture tubes.
- Label one tube “(+)”. Add 10  $\mu$ L of A19 glycerol stock to (+).
- Label the other tube “(-)”. This will be your negative control, testing if your technique is sterile.
- Incubate both tubes overnight shaking at 37°C / (225-250) rpm / (12-16) hr.
- Check if (-) has growth. If not, continue.
- Back dilute overnight 1:250 - 1:1000 into 4x 450 mL fresh LB in 2 L baffled Erlenmeyer flasks (e.g., 1.8 mL overnight into 450 mL LB).
- Incubate back diluted cultures at 37°C / (225-250) rpm to mid-log phase ( $OD_{600}$  between 0.6 and 0.8, typically  $\sim$ 3 hrs).
- Fill 1 L centrifuge bottles with culture. Balance centrifuge bottles against each other and centrifuge cultures at 16 000 rcf / 4°C / 10 min.
- Decant supernatant, add fresh culture, and repeat centrifugation as above, working through the remaining culture. You should end up with large pellets at the bottom of each centrifuge bottle.
- Wash the pellets by resuspending (4°C) NaCl (0.9%) in about 50 mL and transfer the resuspended cells to a single centrifuge bottle. Dilute to  $\sim$ 500 mL and re-pellet at 16 000 rcf / 4°C / 10 min.
- Transfer pellets by spatula into a tared bag weigh and record the mass.
- Flash freeze pellet in liquid nitrogen and store at  $-80^{\circ}\text{C}$ .

### 5. LYSIS

- Resuspend (2-5) g cell pellet in 25 mL of Ribosome Lysis buffer & lyse cells using 130-watt probe sonicator (probe tip diameter: 6 mm) on ice with following parameters: 50% amplitude, 15s on/ 30s off for 2 min on-time. The amount of energy delivered via sonication will vary depending on the amount of cells resuspended.
- Clarify lysate by centrifugation at 16 000 rcf / 4°C / 10 min.
- Aspirate supernatant and measure volume. Add an equal volume of Salting Out buffer to adjust the concentration of ammonium sulfate to 1.5 M and mix well. Incubate at 4°C / 10 min.
- Remove precipitate by centrifugation at 16 000 rcf / 4°C / 10 min.

- Filter supernatant using a 0.22  $\mu\text{m}$  syringe filter and keep cold (4°C).

## 6. FPLC PURIFICATION

- Connect the two Butyl column (5 mL) in tandem, totaling a column volume (CV) of 10 mL.
- Place A1 in Ribosome Wash Buffer and B1 in Ribosome Elution Buffer. Place sample line in A2. Set the default flow rate to 4 mL / min (except for pump washes: 10 mL/ min).
- Perform a pump wash with Ribosome Wash Buffer (without TCEP) and equilibrate the column with 4 CV of Ribosome Wash Buffer(without TCEP).
- Once you've equilibrated your columns, add TCEP to Ribosome Wash and Elution Buffer.
- Load your fraction collector with 15 mL conical tubes and set the fraction volume to 5 mL.
- If using a sample pump to load samples, place sample line (S1) into sample and load around 90% of sample volume onto the column. Once almost loaded, dilute the sample with Ribosome wash buffer (~5 mL) to load as much sample as possible. DO NOT allow air into the FPLC; make sure the sample line is always submerged.
- Wash step 1: wash with 3 CV of Ribosome Wash Buffer to remove unbound components .
- Wash step 2: wash with 5 CV of 80% Wash Buffer and 20% Ribosome Elution Buffer.
- Elution: elute the product by applying 3.5 CV (35 mL) of 50% Ribosome Wash Buffer and 50% Ribosome Elution Buffer. Ensure that the fraction collector captures these fractions separately.
- Wash step 3: Elute all strongly interacting contaminants with 5 CV of 100% Ribosome Elution Buffer.
- Place inlet into a NaOH (0.5 M) and perform pump wash. Wash the column with 3 CV NaOH (0.5M).
- Place the inlet into water, perform pump wash, and then wash column in 2 CV filtered Ultrapure water.
- Place the inlet into AcOH (0.1 M), perform pump wash, and subsequently wash column with 3 CV AcOH (0.1 M).
- Pump wash with water and wash column with 2 CV filtered MilliQ water.
- Place all inlets into EtOH (20% v/v). Perform a pump wash, then wash columns with 3 CV EtOH (20% v/v). Store columns at 4°C in EtOH (20% v/v) until ready for use.

## 7. ULTRACENTRIFUGATION

- Gently overlay recovered fractions (should correspond to second peak) onto 35 mL of Cushion Buffer in a polycarbonate ultracentrifuge bottle (70 mL).
- Prepare another polycarbonate ultracentrifuge bottle as a balance. Measure 35 mL of Cushion Buffer, then add Ribosome Buffer until the balance mass is within 0.1 g of the sample bottle mass. **Make sure all bottles are well balanced ( $\Delta m \leq 0.1$  g) and have no cracks!**
- Pellet ribosomes by ultracentrifugation at 100 000 rcf / 4°C / 16 hrs. A translucent ribosome pellet will be formed at the bottom of the centrifuge bottle. It may be difficult to see.
- Discard the supernatant. Carefully, wash each pellet with 0.5 mL cold ribosome buffer. Repeat this step twice.
- Resuspend the clear pellets in 100  $\mu\text{L}$  of Ribosome Buffer on ice using a magnetic stir bar (3 mm diameter, 10 mm length) on a magnetic stirrer set at the lowest possible speed. Collect resuspended ribosomes.

- Wash tubes with an additional 50  $\mu\text{L}$  of Ribosome Buffer to resuspend any remaining ribosomes.

## 8. QUALITY CONTROL

- Determine the ribosome concentration by measuring the absorbance at 260 nm at a 100x dilution in Ribosome Buffer. 10 units of  $A_{260}$  from a 100x dilution corresponds to 23  $\mu\text{M}$  of undiluted solution.
- Dilute to final stock of 10  $\mu\text{M}$ . To adjust the concentration, dilute the ribosomes with ribosome buffer or concentrate further via centrifugation at 4000 rcf in a 100 kDa centrifugal filter at 4°C.
- Protein gel: dilute 10  $\mu\text{M}$  sample by 4x (Add 2.5  $\mu\text{L}$  of sample with 7.5  $\mu\text{L}$  water) and mix with 2  $\mu\text{L}$  of 6x Laemmli loading buffer. Boil samples at 90°C for 10 min and load 5  $\mu\text{L}$  and 2.5  $\mu\text{L}$  onto 4-20% tris-glycine gel. Run gel at 200 V / 30-45 min or until the loading dye line reaches the bottom of the gel.

## 9. STORAGE

- Aliquot your ribosomes to reduce freeze / thaw cycles and store at -80°C.